

Immunohistochemistry on cryosections using TSA Plus Fluorescence kits

(adapted from *TSA Plus Fluorescence Kit Manual*)

1. Slide Preparation

- a) Air dry slides in oven for 30 min -1 hour.
- b) Wash the slides in 1XPBS for 5 min.
- c) Quench endogenous peroxidase activity by incubating the slides in 3% H₂O₂ in PBS for 10 min.
- d) Wash the slides in 1XPBS for 5 min.

2. Antigen retrieval (Heat-based antigen retrieval)

- a. Preheat the steamer
- b. Preheat the antigen unmasking working solution in a microwave for 10 mins or by heating the solution over a Bunsen burner until boiling
Note: prepare 450ml antigen unmasking working solution (4.2ml antigen unmasking stock solution + 450 ml dH₂O)
- c. Transfer the boiled antigen unmasking solution to the steamer and steam the slides in the antigen unmasking solution for ~15-20mins
- d. After steaming is complete, transfer the slides to FRESH 1XPBS solution

3. Blocking Step

Cover tissue sections with TNB blocking buffer and incubate slides in a humidified chamber for 60 minutes at room temperature.

TNB Blocking Buffer

0.1 M TRIS-HCl, pH 7.5

0.15 M NaCl

0.5% Blocking Reagent

4. Primary Antibody Incubation

- a) Drain off the blocking buffer
- b) Dilute the primary antibody in TNB blocking buffer (concentration depends on the antibodies used)
- c) Add primary antibody to the sides. Use enough volume to completely cover the tissue section (generally 100-300 μ L per slide).
- d) Incubate the primary antibody for 60 minutes at room temperature.

5. Wash

Wash the slides 3X for 5 minutes each in TNT Buffer at room temperature.

TNT Wash Buffer

1M TRIS-HCl, pH 7.5

0.15 M NaCl

0.05% Tween®20

6. Introduction of HRP

- a) Drain off the blocking buffer
- b) Dilute the HRP labeled secondary antibody in TNB blocking buffer at 1:200-1:500
- c) Add primary antibody to the sides. Use enough volume to completely cover the tissue section (generally 100-300 μ L per slide).
- d) Incubate the slices in HRP labeled secondary antibody for 1-2 hours at room temperature.

7. Wash

Wash the slides 3X for 5 minutes each in TNT wash buffer at room temperature

8. Amplification

- a) Pipette 100-300 μL of TSA Plus Working Solution onto each slide.
- b) Incubate the slides at room temperature for 3 to 10 minutes.

9. Wash

Wash the slides 3X for 5 minutes each in TNT wash buffer at room temperature.

10. Visualization of Deposited Fluorophores

For fluorescent detection, counterstain with DAPI and mount for fluorescence microscopy.