

## Immunofluorescence for paraffin-embedded tissue sections

### DAY 1

#### 1) Remove paraffin and rehydrate the tissue sections

- a) Xylene -----2x 5min RT (\*\* Perform this step in the hood\*\*)

Note: Air dry xylene slides in the hood and circle the tissue sections with

Hydrophobic Barrier PAP Pen

- b) 100% Ethanol----- 2-3 min RT
- c) 100% Ethanol----- 2-3 min RT
- d) 100% Ethanol-----2-3 min RT
- e) 95% Ethanol-----2-3 min RT
- f) 70% Ethanol -----2-3 min RT
- g) 50% Ethanol-----2-3 min RT
- h) 1XPBS----- 5 min RT

#### 2) Antigen retrieval (Heat-based antigen retrieval)

- a) Preheat the steamer
- b) Preheat the antigen unmasking working solution in a microwave for 10 mins or by heating the solution over a Bunsen burner until boiling

Note: prepare 450ml antigen unmasking working solution (4.2ml antigen unmasking stock solution + 450 ml dH<sub>2</sub>O)

- c) Transfer the boiled antigen unmasking solution to the steamer and steam the slides in the antigen unmasking solution for ~15-20mins

d) After steaming is complete, transfer the slides to FRESH 1XPBS solution

### 3) Block unspecific staining

- a) Add blocking buffer to the sides (~ 200µl per slide). Make sure the blocking buffer have covered all the specimen
- b) Incubate the tissue sections with blocking buffer (stored in the freezer) inside a humid chamber for 1h at RT
- c) Proceed directly to step 4. **\*\*DO NOT WASH\*\***

### 4) Primary antibody

- a) Dilute the primary antibody in blocking buffer at concentrations depend on the antibodies used
- b) Add primary antibody to the sides (~ 150µl)
- c) Incubate the tissue sections with primary antibody O/N at 4<sup>0</sup>C inside a humid chamber
- d) Wash in PBST 5'X3

Note: Prepare 500ml PBST (50 ml 10xPBS + 0.5 ml Tween 20 + 449.5 ml dH<sub>2</sub>O)

## DAY 2

### 5) Secondary antibody **\*\*\*\*LIGHT SENSITIVE\*\*\*\***

- a) Dilute the secondary antibody in blocking buffer (1:200)
- b) Add secondary antibody to the sides (~ 150µl per slide)
- c) Incubate the tissue sections with the secondary antibody inside a humid chamber for 2h at RT

d) Wash the slides in PBST 5' (COVER WITH ALUMINUM FOIL)

**6) DAPI staining \*\*\*\*LIGHT SENSITIVE\*\*\*\***

a) Add DAPI to the slides (~ 200 $\mu$ l per slide)

b) Incubate the tissue sections with DAPI (stored in the fridge) for 5 mins inside the humid chamber

c) Wash the slides in PBST for 5min (COVER WITH ALUMINUM FOIL)

d) Mount the slides for imaging