

Whole mount immunofluorescence on embryonic tissue slices

(adapted from Gómez-Gaviro et al., 2017)

DAY 1

1) Fix embryonic tissue slices (~200-400µm in thickness) in formalin

- a. Incubate the samples in 10% formalin O/N at RT
- b. Wash in PBST 5'X3

Note: Prepare 500ml PBST (50 ml 10xPBS + 0.5 ml Tween 20 + 449.5 ml dH₂O)

DAY 2

2) Reagent 1 ******LIGHT SENSITIVE******

- a. Prepare Reagent 1 in 50 ml tube, and shake until dissolved:
 - i. Distilled water 14 ml
 - ii. Urea 10 g
 - iii. N, N, N, N – tetrakis ethylenediamine 10 ml
 - iv. Triton 6 ml
 - v. DAPI 1:5000
- b. Incubate the samples in Reagent 1 O/N at 37⁰C, rotating. Make sure the samples are fully immersed in Reagent 1 (depends on the size of the sample)
- c. Wash in PBST 5'X3

3) Antigen retrieval (Heat-based antigen retrieval) ******LIGHT SENSITIVE******

- a. Preheat the steamer

- b. Preheat the antigen unmasking working solution in a microwave for 10 mins or by heating the solution over a Bunsen burner until boiling

Note: prepare 450ml antigen unmasking working solution (4.2 ml antigen unmasking stock solution + 450 ml dH₂O)

- c. Transfer the boiled antigen unmasking solution to the steamer and steam the samples in the antigen unmasking solution for ~10-15mins
- d. After steaming is complete, transfer the samples to FRESH 1XPBS solution

4) Block unspecific staining ******LIGHT SENSITIVE******

- a. Incubate samples in blocking buffer for 4 hours at 37⁰C, rotating. Make sure the samples are fully immersed in blocking buffer (depends on the size of the sample)
- b. Proceed directly to step 5. ****DO NOT WASH****

5) Primary antibody ******LIGHT SENSITIVE******

- a. Dilute the primary antibody in blocking buffer (concentration depends on the antibody used)
- b. Incubate the samples with primary antibody 4 hours at 37⁰C, rotating. Make sure the samples are fully immersed in primary antibody (depends on the size of the samples)
- c. Wash in PBST 5'X3

Note: Prepare 500ml PBST (50 ml 10xPBS + 0.5 ml Tween20 + 449.5 ml dH₂O)

6) Secondary antibody ******LIGHT SENSITIVE******

- a. Dilute the secondary antibody in blocking buffer (1:200)

- b. Incubate the samples with the secondary antibody for 4 hours at 37⁰C, rotating.
Make sure the samples are fully immersed in secondary antibody (depends on the size of the sample)
- c. Wash the slides in PBST 5'X3

7) Reagent 2

- a. Prepare Reagent 2, shake in 50 ml tube until dissolved:
 - i. Distilled water 1.5 ml
 - ii. Urea 2.5 g
 - iii. Sucrose 5 g
 - iv. Triethanolamine 1 ml
- b. Incubate the samples in Reagent2 O/N at 37⁰C, rotating
- c. Mount the samples for imaging

Reference:

Gómez-Gavero, M.V., Balaban, E., Bocancea, D., Lorrio, M.T., Pompeiano, M., Desco, M., Ripoll, J., and Vaquero, J.J. (2017). Optimized CUBIC protocol for three-dimensional imaging of chicken embryos at single-cell resolution. *Development* 144, 2092-2097.